

Gender Detection System - User Guide

Version: 1.0

Last Updated: January 2025

Table of Contents

1. [Introduction](#)
 2. [System Requirements](#)
 3. [Installation & Setup](#)
 4. [Using the System](#)
 - [Uploading Images](#)
 - [Viewing Detection Results](#)
 - [Generating Reports](#)
 - [Managing Reminders](#)
 - [Updating Settings](#)
 5. [Troubleshooting](#)
 6. [Frequently Asked Questions \(FAQs\)](#)
 7. [Support](#)
-

1. Introduction

The Gender Detection System is a web-based application designed to analyze **teeth radiographs** and determine the gender of individuals using an AI-powered **YOLOv11 model**. It provides insights useful for:

- **Dentists**
- **Forensic Experts**
- **Medical Diagnosticians**

The system allows users to: ☒ Upload radiograph images (single or batch).

☒ View gender detection results.

☒ Generate visual reports with confidence scores.

☒ Set reminders for follow-ups or reviews.

2. System Requirements

Hardware:

- Minimum **4GB RAM** (8GB+ recommended)
- **CPU:** Intel i5 or higher (**GPU acceleration optional**)
- **Storage:** At least **5GB free space**

Software:

- **Python 3.8+**
 - **Flask Web Framework**
 - **Ultralytics YOLOv11**
 - **OpenCV**
 - **Matplotlib**
 - **FPDF for Reports**
-

3. Installation & Setup

Step 1: Install Dependencies

Run the following command to install required libraries:

```
bash
CopyEdit
pip install flask ultralytics opencv-python matplotlib fpdf werkzeug
```

Step 2: Download and Setup Model

Ensure your **trained model file (best.pt)** is placed inside the `models` folder:

```
bash
CopyEdit
/models/best.pt
```

Step 3: Run the Application

Navigate to the project directory and run:

```
bash
CopyEdit
python app.py
```

The system will start and be accessible at:

4. Using the System

Uploading Images

1. Click on the **Upload Data** section.
 2. Choose one or multiple radiograph images (JPG, PNG).
 3. Click **Upload**.
 4. The system will process the images and display:
 - Detected **gender**.
 - **Confidence score**.
 - **Detection visualization**.
-

Viewing Detection Results

Once the images are processed:

- A **summary** of male vs. female detections will be displayed.
 - **Low-confidence detections** will be flagged.
 - A **graph visualization** of results will be generated.
-

Generating Reports

1. After uploading images, navigate to **Reports**.
 2. Apply filters:
 - Filter by **Male / Female**.
 - Set a **Minimum Confidence Score**.
 3. Click **Download Report**.
 4. A **PDF report** will be generated with:
 - Detection results.
 - Confidence scores.
 - A visual summary chart.
-



Managing Reminders

1. Go to **Reminders**.
 2. Enter:
 - **Task Description**
 - **Date & Time**
 3. Click **Add Reminder**.
 4. Active reminders will be displayed with a delete option.
-



Updating Settings

1. Navigate to **Settings**.
 2. Update:
 - **Email**
 - **Password**
 3. Click **Save Changes**.
-

5. Troubleshooting

Issue	Solution
Images not uploading	Ensure images are in JPG/PNG format.
Model not detecting gender	Ensure <code>best.pt</code> is in the correct location.
Report not generating	Check if results exist before clicking Generate Report .

6. Frequently Asked Questions (FAQs)

Q1: Can I upload multiple images at once?

Yes, the system supports **batch processing** of multiple images.

Q2: Can I filter detection results?

Yes, in the **Reports** section, you can filter results based on **gender and confidence scores**.

Q3: What if the confidence is low?

The system **flags low-confidence detections**, allowing for **manual review**.
